

# » TECHNICAL DATA

	Club Sport	Elite	Grand Slam	TIGA	TIGA fx	HiLite
<b>Seat width:</b>	340, 370, 400, 440 mm	250 to 500 mm	250 to 500 mm	250 to 500 mm	250 to 460 mm	250 to 500 mm
<b>Seat depth:</b>	300, 330, 360, 390, 430 mm	250 to 500 mm	250 to 500 mm	250 to 500 mm	250 to 460 mm	250 to 500 mm
<b>Front seat height:</b>	480, 500, 530 mm	360 to 630 mm	360 to 630 mm	380 to 550 mm	480 to 550 mm	380 to 550 mm
<b>Rear seat height:</b>	400, 420, 440, 530 mm	360 to 630 mm	360 to 630 mm	360 to 550 mm	360 to 550 mm	360 to 550 mm
<b>Backrest height:</b>	Adjustable from 250 to 300 mm	150 to 450 mm	150 to 450 mm	150 to 600 mm	150 to 400 mm	150 to 600 mm
<b>Backrest angle:</b>	Fixed at 0°	From -10 to +10	De -10 to +10°	From -10 to +20	From -7 to +7	From -10 to +20
<b>Camber:</b>	18 °	16, 18, 20, 22 °	16, 18, 20, 22 °	0, 1, 2, 3, 5 °	0, 1, 2, 3, 5 °	0, 1, 2, 3, 5 °
<b>Max. user weight:</b>	125 kg	125 kg	125 kg	125 kg	125 kg	125 kg
<b>Product weight:</b>	10 to 13.5 kg	9.5 to 12.5 kg	9.5 to 12.5 kg	6.5 to 11.5 kg	8.5 to 11.5 kg	7.5 to 11.5 kg
<b>Frame colours:</b>	23	25	25	25	25	25

For further information on the full specification, options and accessories, please refer to the order form.  
All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



## Sunrise Medical GmbH & Co. KG

Kahlbachring 2-4  
69254 Malsch / HD  
Germany  
Tel: +(49) (0) 7253 980 0  
Fax: +(49) (0) 7253 980 222  
[international@sunrisemedical.de](mailto:international@sunrisemedical.de)



**BELIEVE | LIVE | INSPIRE**

# » RGK As individual as you are

More than three decades ago the RGK founders meet through their mutual passion for wheelchair basketball.

Together we identified using the same chair for general life as our favourite sport was not “cutting the mustard”, so we investigated rigid frame chairs that had newly arrived to market and discovered there was still nothing of a truly bespoke, tailored design available. The idea for RGK was born

Thirty years on and we have an amazing team. We are thrilled and excited that their skill, passion and commitment have turned our idea into a brand with a global presence and a reputation for excellence.

# A Jama

Abdi Jama, GB men's wheelchair basketball captain  
2x Paralympian Bronze Medallist  
2x European GOLD Medallist

Respect each person individuality is our RGK philosophy. Because we believe that every person is unique. That is why our chairs are created and built only for you. Fully bespoke to adapt to your measures, body shape, your preferences, your life.

This complete customization allows us to offer you the best quality chair, that reacts immediately to each of your movements. Discover with RGK a world of sensations. The perfect transmission of your energy.



# » AS INDIVIDUAL AS YOU ARE

A chair that follows your body shape and adapts to your precise measures, like a second skin.

A light chair that reacts to each of your movements. A world full of sensations. a direct transfer of your energy!

With maximum design and quality, and tailored to the minimum detail to reflect your personality.

The perfect combination.



# HOW YOUR CHAIR IS CREATED



## 1. INTRODUCTION

Understanding you, talking through your needs and lifestyle and taking your precise measurements gives us the benchmark to initiate our design process.



## 2. TAKING YOUR MEASUREMENTS

We take your precise measurements to make sure the chair will perfectly follow your body shape and fit you like a second skin.



## 3. DESIGN

Your personal measurements are introduced to our CAD system giving accurate 3D imagery, paving the way for a flawless creation process.



## 4. CREATION

Your chair is brought to life, handmade by our experienced team of technicians.



## 5. BUILD

With the component parts completed, assembly, rigorous testing and inspection are carried out.



## 6. YOUR PERSONALIZATIONS

Your chair is an extension of you so you can tailor the look of your chair to your personal taste, preference and personality from frame colours to any kind of upholstery embroidery.

# TIGAF~~X~~

» The rigid chair that folds



# For people on the move!

The custom chair that meets perfectly your daily and travel needs.

Airline travel? Small car?, No problem! TIGA FX folds into a briefcase style shape, ideal for storage in tight spaces. The TIGA FX easily fits into the overhead luggage compartment and the smallest car boots! Our Q-lock system ensures all the benefits of a folding chair without compromising strength. The system is simple, clean and streamlined enough to meet all your travel needs. Its frame is created individually for you, fully customized to your measures, needs and personal preferences.

## Key Features


- » Aluminium 7020 customized frame
- » Q Lock folding frame and backrest
- » Integrated fix COG or Adjustable COG available
- » Ergonomic seat and back optional
- » Fold Flat sideguards
- » Integrated footrest with carbon fibre cover option
- » Wide accessories range
- » Choice of anodised accent colours
- » Airtech upholstery system
- » Max user weight 125 kg



Q-lock system for an easy folding of front frame and backrest.



"Fold flat" folding sideguards.



Ergonomic seat and back  
Optimizes user posture to increase stability, facilitate propulsion, and help skin protection.



## FRONT WHEEL

Enables you to travel easily across all sort of terrains by keeping your front castor wheels off the floor at all times. Attaches to TIGA, TIGA FX and Hi Lite.





Be ready to discover new worlds with the Frontwheel!

Add it to your RGK chair in seconds and enjoy the freedom of movement in any kind of terrains.

Great adventures are waiting for you. Go for them!

# Pushing the boundaries, see where it will take you

The frontwheel enables you to travel across grass, cobbles, mud, gravel and other uneven ground easily by keeping your front castor wheels off the floor at all times when you move in any direction, including backwards.

## Technical Specifications

- » Exclusive to the TIGA and TIGA fx
- » Three point frame connection
- » Compact and easy to install in seconds
- » Fits to frame with support under footrest for strength
- » Castor wheels remain off the ground at all times in any direction
- » Footrest width is unaffected
- » No need to move feet during installation
- » Locks up and down for when both in use and not in use



200x50 pneumatic tyre with anti-flutter fork.



Quick release handles to activate and release.



Robust three point connection to frame.



Retro fit to any TIGA/TIGA fx via mounting clamps.