

SmartDrive MX2+

Independence rediscovered



permobil

Power every adventure

SmartDrive's intuitive controls add a power assist for a variety of manual wheelchair scenarios like cruising up hills or over thick carpets. Its lightweight design attaches quickly and its unique omni-wheel design moves with you as you pivot, turn, wheelie, and navigate curbs in style.



*Extensive battery life for all day usage on a single charge

PushTracker

Use the PushTracker band and Smartphone app to effortlessly control your SmartDrive and track your daily activity. By tracking your push count, coast time, and daily distances, the PushTracker app can help you monitor your efficiency, give you a more complete picture of your activity, and make individualized adjustments.



Compatible with iOS 10.0 / Android 5.0 and newer.





Easy on, easy off

The SmartDrive is lightweight and easy to install on virtually any make and model of manual wheelchair including folding, tilt-in space, one-arm drive and standing. At only 12.5 lbs (5.7 kg), the MX2+ is easily attached or removed from the back of your chair, making it the perfect on-demand addition to your mobility.

Spring loaded flanges for quick attachment on folding chairs.

Specifications

Product Dimensions	Length: 15.3 in (389 mm) Width: 5.7 in (141 mm) Height: 9.5 in (242 mm)
Product Weight	12.5 lbs (5.7 kg) Drive Unit
User Weight	30 to 331 lbs (14 to 150 kg)
Driving Range*	12.3 miles (19.8 km) Level/No Resistance
Motor Power	250 W Brushless DC
Maximum Speed	5.5 mph (8.8 km/h)
Operating Temp	-13°F to 122°F (-25°C to 50°C)
Chair Type Compatibility	Rigid, Folding, Tilt-In Space, One-Arm Drive and Standing
Drive Wheel Diameter	22, 24, 25, 26 in (501, 540, 559, 590 mm)
Warranty	2-Year Limited Warranty
Bluetooth® LE	FCC ID QOQBT113 IC ID 5123A-BGTBLE113
HCPCS	E0986

**Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.

A wheelchair's specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.



Please refer to your owner's manual for more detailed operation and care instructions.

Since U.S. Medicare coding and coverage criteria is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.
©2018 Permobil. Permobil has a policy of continual product improvement and reserves the right to amend this document.